

Mycoplasma pneumoniae (MP)-IgM Antibody

(Technology of colloidal gold immunochromatography)



QUICK DETAIL

Sensitivity: 99.2 %

Specificity: 99.5 %

Read Time: 15-20 minutes

Method: Immunochromatography test
(colloidal gold)

Specimen: serum or plasma

Certificate: CE Marking

SPECIFICATIONS

1. For In vitro diagnostic use only
2. Used as quantitative screening test
3. Package 20 Tests/Kit



FEATURES

- ✓ Easy to use
- ✓ Easy to read
- ✓ High accuracy
- ✓ High speed